

# Prism SA-200

## Aluminium Stand Aid



## Owner's Manual

Use and Care  
Trouble Shooting  
Warranty Information

Introduction.....	2
Overview.....	2
Components of the Stand Aid.....	3
Component List .....	4
Specifications .....	5
Cautions .....	6
Assembly Instruction .....	7
Operation	
Charging.....	12
Start.....	12
Forward movement.....	13
Turning.....	13
Leg Adjustment.....	14
Raising and Lowering Boom.....	14
Castors and Braking.....	14
Basics in transferring an individual.....	15
Emergency Stop & Lowering.....	17
Control Box-Usage Information.....	18
Control Box-Service Information.....	19
Technical Specifications .....	20
Do's and Don'ts .....	21
Fault finding.....	22
General inspection and maintenance.....	23
Points of Attachment - Checklist .....	24
Service record history.....	25
Warranty.....	29

**PRISM MEDICAL UK**

Prism Medical UK  
Unit 4 Jubilee Business Park  
Jubilee Way  
Grange Moor  
Wakefield  
WF4 4TD

Tel: 0845 409 2000  
Fax: 0844 980 2261  
Website: [www.prismmedical.co.uk](http://www.prismmedical.co.uk)



**CAUTION: DO NOT ATTEMPT TO USE THIS EQUIPMENT WITHOUT FIRST UNDERSTANDING THE CONTENTS OF THIS MANUAL.**

## Introduction

Before using this equipment, and to ensure the safe operation of your Prism SA-200 Stand Aid, carefully read this entire manual, especially the section on “Cautions”. The SA-200 is designed to be used in conjunction with Prism Medical UK accessories and slings. Please refer to any user guides supplied with these components and refer to them while reviewing this manual.

Should any questions arise from reviewing this manual contact your local authorised Prism Medical UK representative. Failure to comply with warnings in this manual may result in injury to either the operator, or the individual being lifted/transferred. Damage to the Stand Aid and/or related components may also occur. Be sure that the contents of this manual are completely understood prior to using this Stand Aid.

Store this manual with the documents included with the stand aid and sling (s). Contents of this manual are subject to change without prior written notice.

---

## Overview of Prism SA-200 Stand Aid

The Prism SA-200 Stand Aid is a lifting aid used by health care professionals to transfer clients. The Stand Aid makes it possible to move mobility impaired individuals with minimal strain or risk to the caregiver, while providing complete safety, dignity and comfort for the client.


The Prism SA-200 Stand Aid is one of two components that make up this technology. The other component, the sling, is a specially designed fabric accessory that attaches to the stand aid by means of a sling hook and straps, and holds an individual while the lift, or transfer takes place. The sling is generally supplied with the Stand Aid at the initial time of purchase. Please refer to any user guides supplied with the sling and reference them while reviewing this manual.

The Prism SA-200 Stand Aid has the ability to stand an individual up from one location such as chair, then move the individual to another location and finally lower the individual into a wheelchair. The functions of lifting up or down are accomplished by pressing buttons on the hand control. The hand control is attached to the Stand Aid. Due to the design of the Stand Aid system, it takes very little effort to press a button to perform the desired motion.

Please familiarise yourself with the components of the Prism SA-200 Stand Aid by referring to the diagram on the next page.

<b>Model Table for Prism SA-200 Stand Aid</b>	
<b>Product Code</b>	<b>Product Description</b>
SA-200	SA-200 - Aluminum Stand Aid 200 Kgs

# Components of the Prism SA-200 Stand Aid



# Component List

The following components are included with your new Prism SA-200 Stand Aid system:

- Prism SA-200 Stand Aid
- Hand Control
- Stand Aid Integrated Charger
- Mains Charger Cable
- Owner's Manual
- Multi Function Spanner



- 5 mm Hex Allen Key (Short Arm)



- 4 mm Hex Allen Key (Short Arm)



- 3 mm Hex Allen Key (Short Arm)



**SLINGS:** If a sling has been supplied with the Stand Aid, then please refer to the instructions included with the sling.

**IMPORTANT:** Before initial use, the Stand Aid unit must be charged for 4 to 5 hours. Refer to section titled “**Charging Instructions**”. The Hand control must also be connected to the Stand Aid. If it is not connected, then please refer to the section titled “**Assembly Instructions**”.

# Technical Specifications- SA-200 Stand Aid

## (A) Dimensions and Weight - SA-200 Stand Aid

Specification	SA-200
Maximum weight capacity	200 kg (31.5 Stone)
Maximum lifting height	138 cm
Minimum lifting height	90 cm
Height to top of legs	11 cm
Clearance from bottom of legs to floor	5 cm
Overall length	103.5 cm
Distance inside the legs (min.)	58 cm
Distance inside the legs (max.)	89 cm
Distance outside the legs (min.)	68 cm
Distance outside the legs (max.)	99 cm
Turning radius	117 cm
Wheels (dual-wheel castors)	Front –76 mm
	Rear – 76 mm with brake
Weight of Stand Aid	36.8 kg
Total Shipping Weight with Box	43.83 kg

## (B) Electrical Specifications - SA-200 Stand Aid

<b>Lift Motor:</b>	24 VDC, 10.5 Amps Max.
<b>Charger Input:</b>	100-240 VAC, 50-60 Hz, 0.4 Amps Max, Class 2, IPX4
<b>Charger Output:</b>	24 VDC, 0.65 Amps
<b>Batteries:</b>	24 VDC (2 x 12 VDC) 2.9 Ah, Sealed Lead Acid
<b>Hand Control:</b>	Electric, IPX4
<b>Lifting Range:</b>	900mm to 138mm
<b>Stand Aid Weight:</b>	36.8 Kg
<b>Maximum Load:</b>	200 Kg (31.5 stone)
<b>Duty Cycle:</b>	10% use, 90% rest
<b>Rated Performance:</b>	30 lifts at 440 lbs., 10% duty cycle.

**Maximum load of the installed Stand Aid is determined by referring to the product label located on side of the Stand Aid**

### Shipping/Storage Conditions:

Temperature: -40 to +70 °C  
 Relative Humidity: 10 to 100% RH  
 Atmospheric Pressure: 500 to 1060 hPa

### Normal Operating Conditions:

Temperature: +10 to +70 °C  
 Relative Humidity: 30 to 75% RH  
 Atmospheric Pressure: 700 to 1060 hPa



## Cautions:

- Under no circumstance should the Prism SA-200 Stand Aid and slings be put in control of a person who has not been properly trained in the use and care of this equipment. Failure to adhere to this warning may result in serious injury to the operator, and/or the individual being lifted/ transferred.
- The Prism SA-200 Stand Aid and slings are not toys. Do not use it for unsafe practices. Do not allow children to play with the Stand Aid or any of its components.
- The manufacturer's warranty is voided if persons unauthorised by Prism Medical UK perform work on the Prism SA-200 Stand Aid.
- There are no user serviceable parts inside the actuator. Do not remove cover screws, or open the unit, as this will **VOID THE WARRANTY**.
- In facilities where more than one operator will be responsible for using the Prism SA-200 Stand Aid and sling (s), it is imperative that all such members are to be trained in its proper use. A training programme should be established by the facility to acquaint new operators with this equipment.
- Never expose the Prism SA-200 Stand Aid directly to water. Warranty does not cover any misuse or abuse of the Stand Aid system.
- To maintain optimum function, the Prism SA-200 Stand Aid should be inspected and maintained on a regular basis. See the section titled “General Inspection and Maintenance”.
- Any accessories used with the Prism SA-200 Stand Aid including sling (s), should be checked to ensure that they are in good working order. Check for signs of wear or fraying prior to use. Report any unusual wear, or damage immediately to your local authorised Prism Medical UK Service Provider.
- The Prism SA-200 Stand Aid and associated slings are intended **only** for lifting and transferring of a person. Prism Medical UK will not be responsible for any damage caused by the misuse, neglect or purposeful destruction of the Stand Aid, and/or its associated components.
- In any circumstances do not exceed the maximum allowable load of this Stand Aid. Refer to the “Specifications” section of this manual, and/or the labels on the side of the Stand Aid.
- There is a risk of explosion if the Stand Aid is used in the presence of flammable anaesthetics.
- Ensure that a clear space is maintained around the Stand Aid. Move any obstacles out of the way before operating the Stand Aid.

## Assembly Instructions:

1. Place the Carton/Packaging box in a clear working area, open carefully, and remove following items and place on the floor, taking care to protect the finish from damage.
  - Mast with Boom assembly
  - Base Assembly with Legs attached
  - Battery Pack, Control Box and Charging Power Cord
  - Hand Control
  - Tool Kit
  - User Manual, and Test Certificate



**Caution: Component parts are heavy, and will need to be lifted with care.  
Heavier items may need two people for assembly.**

2. Place the Base assembly with legs attached on the floor, and lock both rear castors as shown in Figure 1.



**Figure 1**

3. Mount Knee Pad Height Adjustment Rod using four screws supplied with the Stand Aid as shown in Figure 2, and 3.



**Figure 2**



**Figure 3**

## Assembly Instructions:

- Place the Mast with Boom assembly upright in the tube on the Base assembly as shown in Figure 4. Using Locking Handles, secure the mast in the desired position on the base as shown in Figure 5. Adjust the position of the locking handles pointing in downward direction towards the base as shown in Figure 4.



Figure 4




Figure 5



**Caution: Possible finger trap. Keep fingers away from the end of the mast when inserting into the tube on the Base assembly. Heavier items may need two people for assembly.**

- Route the Lifting Actuator Cable assembly as shown in the Figure 6.



Figure 6



## Assembly Instructions:

6. Pull the plastic knob to insert the Knee Pad Assembly into the Knee Pad Height Adjustment Rod and release the plastic knob once it aligns with hole on Knee Pad Height Adjustment Rod as shown in Figure 7, and 8.



**Figure 7**



**Figure 8**

7. Insert the Sling Bar into the Boom of the Stand Aid by pulling plastic knob on top of the boom as shown in Figure 9. Make sure that the sling bar locking screw is out from the boom as shown in Figure 10. After inserting the sling bar, tighten the sling bar locking screw so that sling bar secured into the boom.



**Figure 9**



**Figure 10**



**Caution: Sling bar locking screw must be tight, and locked for all the times.  
Failure to do so, may result in severe injury.**

## Assembly Instructions:

8. Remove the screws from the mounting position of push handle on Mast assembly as shown in Figure 11. Insert the screws through the push handle bracket, and mount the push handle on the Mast Assembly as shown in figure 12.



**Figure 11**



**Figure 12**

9. Remove the middle screw from the mounting bracket using 3 mm Allen Key. Press the SA-200 CBJ Care Control Box against the mounting bracket and slide it down. The control box should be seated securely on top of the edge of bracket as shown in Figure 13. Mount the Control Box to the mounting bracket using same screw that you removed from the middle of the bracket as shown in Figure 14.



**Figure 13**



**Figure 14**

## Assembly Instructions:

10. To install the battery, you will need to press the battery against the bracket and slide it down until it sits straight on top of the control box as shown in Figure 15. Make sure to secure the battery box by latching the battery lock into the bracket as shown in Figure 16.



Figure 15



Figure 16

11. Connect all the cables to the control box as shown in the Figure 17 to 19.
- A. Connect the Hand Control connector to the large circle on the left side of the control box as shown in the Figure 17.
  - B. Connect the Up/down Actuator Cable connector to the second circle (Marked No. 1) from the left side of the control box as shown in the Figure 18.
  - C. Connect the Mains Power Cord to the rectangular outlet on the control box as shown in the Figure 19.



Figure 17



Figure 18



Figure 19

12. **Please make sure to charge battery before initial use and before each operation of the Stand Aid.**

### **After assembly, check to ensure that:**

1. The mast is fully locked into position.
2. The RED emergency stop button on control box is in the OFF (out) position.
3. Lifting actuator is moving up and down with the buttons on the hand control, and control box.
4. Emergency lowering (mechanical, and electrical) function works properly.
5. Rear wheel brakes works properly.
6. Leg spreading foot pedal works properly.
7. Battery is fully charged.

## Operation:

**⚠ Caution: Visually inspect the Stand Aid before using for any unusual wear and tear. Should anything look unusual then contact your local representative prior to use.**

**Failure to comply with this caution could result in serious injury to the operator, the individual being lifted and/or damage to the Stand Aid.**


**Because of the safety function in the actuator the lifting arm and hanger bar can “fall down” to the height of the spindle. Therefore, both before and after each lifting operation, press the lifting arm down to ensure that it is supported by the actuator. This must be done, to remove the possibility of the lifting arm “falling down” and the risk of subsequent injury occurring.**

## Charging:

The “On” light illuminates when you plug the control box to the power outlet and should turn off when it is unplugged. The “Charge” light illuminates, and when you connect the battery to the control box and it is charging, also the LCD displays shows the symbol of mains cable plug in as shown in the Figure 20. The “Charge” light will turn off when the battery is fully charged.

The Stand Aid should be plugged in for charging whenever it is not in use to get the max number of cycles out of the battery. **An audible battery alarm will sound when approximate 25% of battery capacity remaining.** To obtain max life, charge the battery at no more than 30% discharge.

Please Note: The battery can be charged whether the Emergency Stop is activated or not.



## Start:

The Hand Control has two functions up and down. Press down on the symbols to operate the desired function. It is not possible to use two functions at the same time.

An audible alarm will sound once the load becomes too great on the Stand Aid. This will cause the lifting operation to stop. Do not continue to operate the handset by repeatedly pressing the buttons. At this point the Stand Aid has reached its maximum lifting capacity. This is a purposefully built in safety feature.

## **Please Note!**

**The emergency stop must be released at all times during normal operation.**

# Operation:

## Forward Movement:

To move the Stand Aid forward, hold onto the handle bar and push forward. Alternatively, the caregiver can also place one foot on the base tube and push forward with their foot.



## Turning:

When turning the Stand Aid, stand on along one side. With one hand pull gently on the handle bar, and with the other hand push gently away on the lifting arm. In this way the Stand Aid will rotate around its own axis. This movement is performed with a smooth, slow action to avoid swinging the patient unnecessarily.



## Warning:


Do not hold on to the actuator when turning the Stand Aid as there is a risk of damaging the actuator.

# Operation:

## Leg Adjustment:

The legs of SA-200 Stand Aid are mechanically adjustable for opening, and closing base width. The legs can be opened to enable access around arm chairs or wheel chairs. For transferring through narrow doorways, and passages the Stand Aid legs should be in closed position.

Leg adjustment is achieved by pressing the foot pedal located at the rear of base. To open the legs, press down the left side, and to close the legs press down the right side as shown below.



## Raising and Lowering the Boom:

The up and down movement of boom on SA-200 Stand Aid is achieved by a powerful electric actuator which is controlled by hand control. The hand control has two buttons with directional arrows UP and DOWN. The actuator stops automatically at the limit of travel in both directions.

## Castors and Braking:

The SA-200 Stand Aid has two rear castors with a brake. The rear castors can be braked for rotation, lateral movement, and parking. To apply the brake, press the brake pedal down with your foot. To release the brake, press the raised pedal towards the wheel.

During the lifting, the rear wheels should remain unlocked so that the Stand Aid will move to the patient's centre of gravity. The wheels should, however, be locked if there is a risk of the Stand Aid moving to the patient. For example, when lifting the patient from the floor.



**Warning:** Locked wheels during lifting will increase the risk of the Stand Aid tilting over.

## Application:

**If the Stand Aid is used incorrectly, any warranty or product liability will cease to be valid.**

The Stand Aid must only be used for person lifting and only for persons who, including the sling, do not weigh more than the stated Safe Working Load. If the load exceeds the stated weight limit and if the Stand Aid is used for lifting anything but persons, then any product liability that Prism Medical UK might have, in connection with insurance / warranty / maintenance etc, will cease to be valid.

To avoid possible accidents and injury to persons being lifted, the Stand Aid must only be operated as described in the preceding pages.

## Basics in transferring an individual

**⚠ Caution:** The following steps are intended to generally illustrate the procedure involved in the transferring of an individual from one location to another using the Stand Aid.

The manual for the sling that was purchased with the Stand Aid should be reviewed in detail prior to attempting these steps, as the sling illustrated here may not be the same as the one that was purchased.

Contact your local authorised Prism Medical UK Representatives if you have any questions or concerns.

**Step 1)** Unplug the Stand Aid from its charging station or current location and move close to the individual that is to be transferred. Use the procedures for up and down and transferring as described in the sections titled, “Start”, “Forward Movement” and “Turning”.

**⚠ Caution:** Always use extreme care when moving the Stand Aid from one location to another. Watch out for and avoid any obstructions that may cause injury to the individual in the sling, or damage to the Stand Aid.

**Step 2)** Prepare the individual being transferred with the appropriate sling. Refer to the instructions supplied with the sling that was purchased on how to properly fit an individual with a sling.

**⚠ Caution:** Always make sure that the sling is correctly fitted and adjusted on each side of the individual so that maximum comfort and safety are achieved prior to lifting.

**Step 3)** Once the individual has been fitted with the sling, move the Stand Aid so that it is positioned **directly in front** the individual and utilise the leg opening function if required. Lower the boom to a height so that the straps of the sling can be easily attached to the sling hooks.

**⚠ Caution:** Always check to ensure that the Stand Aid is correctly positioned directly above the person to be lifted.

**⚠ Caution:** Check to ensure that the sling hooks have no cuts, dents or sharp edges that may come in contact with the straps of the sling and cause damage to them. Report any concerns to your local authorised representative.

**Step 4)** Attach the straps of the sling to the sling hooks. Be sure to double check to ensure that the straps are properly attached to the sling hooks, and that the individual being lifted is properly positioned in the sling prior to lifting.

**⚠ Caution:** Prior to lifting an individual make sure that the straps of the sling are securely placed on the sling hooks.

## Basics in transferring an individual ... continued

**Step 5)** The individual may now be raised using the UP button on the hand control. While lifting is in progress the height required in order for the transfer to be completed safely should be closely observed. Ensure that the individual being lifted will not be injured by any obstructions during the initial lifting.

**⚠ Caution: Always use caution when lowering/raising an individual who is in the sling of the Stand Aid. Watch out for and avoid any obstructions that may cause injury to the individual, or damage to the Stand Aid.**

**Step 6)** Once at the correct height the individual can be moved to the desired location. Refer to the sections already described in this manual on how to do so.

**Step 7)** Once at the desired location the individual in the sling can be lowered/raised to the correct height in order to complete the transfer. On completion of lowering/raising ensure that the individual is properly positioned and safely supported prior to removing the straps of the sling from the carry bar.

**⚠ Caution: Prior to removing the straps of the sling from the carry bar be sure to check that the individual being lifted is securely supported in the final desired position.**

**Step 8)** Lower the carry bar sufficiently to allow the straps of the sling to be easily removed from the carry bar. Take care not to let the carry bar come in contact with the individual in the sling. The straps from the sling can now be removed from the carry bar. The Stand Aid should then be moved away from the immediate area so that it will not interfere with the removal of the sling from the client.

**Step 9)** The sling can now be gently removed from the individual. It should then be stored in a safe place for future use.

**Step 10)** The Stand Aid can now be moved to a safe location until further use. The Stand Aid should be turned off when not in use. It is recommended that the Stand Aid be left on charge when not in operation. Refer to the section titled, "Charging the Stand Aid" for instructions on charging.



## Emergency Stop & Lowering:

**⚠ Caution:** The manual emergency lowering system should be used only if the lowering procedures described in the previous section of the manual do not work. Should you have any concerns or questions, contact your local authorised Prism Medical UK Representative.


**⚠ Caution:** DO NOT use the Stand Aid after the manual lowering mechanism has been used. The Stand Aid must be reset by a qualified lift technician after use. Contact your local authorised Prism Medical UK Service Provider.

### Emergency Stop:

Push the Emergency Stop button to cut all power on the Stand Aid. To resume power, release the emergency stop button by turning clockwise as indicated by arrows on the Emergency Stop button.

### Electrical emergency lowering/raising:

The Stand Aid boom is lowered or raised by pressing the UP and DOWN arrow button on the control box.



### Mechanical emergency lowering:

In case of power failure, it is possible to mechanically lower a patient placed in SA-200 Stand Aid. In order to lower the boom in an emergency situation, pull up the red lever located directly on the Stand Aid's actuator. The boom will lower as you press up and hold the lever. Release the red lever once you have lowered the boom to a safe position.



## Control Box: Usage information

The control box of **SA-200** Stand Aid has an LCD display showing basic usage information on the display. The following information will be visible on the LCD display of control box.



When it is time for service the service symbol will appear on the display will light up and give notice to users about the need for service.



If the lift stops because of overload, the symbol will appear on the display and gives the user the information about why the Stand Aid has stopped.



If someone tries to use the Stand Aid whilst plugged in, the symbol will appear on the display. The user receives the information about why the Stand Aid is not working. The use of the Stand Aid whilst charging is not possible.



The display showing a full battery symbol indicates that the battery is fully charged.



The display showing a half empty battery symbol indicates that it is time to charge the battery.



The display showing a empty battery symbol indicates that the battery is completely discharged.



The display showing a lifting actuator is moving up.




The display showing a lifting actuator is moving down.

## Control Box: Service information

The control box of SA-200 Stand Aid has an LCD display showing basic usage information on the display.

To get this information, press the “lifting arm up” button on hand control for half a second.

The information that appears on the display is as shown below.



### **Total cycles made by lifting actuator:**

One cycle is defined as; Moving with load (the actuator draws more than 1.5 Amps).

Moving direction up for a minimum of 5 seconds (several activations are allowed), followed by moving down for a minimum of 2 seconds.

### **Total work made by lifting actuator:**

Work indicator for the lifting actuator measures via ampere usage x seconds in use.

The work indicator gives very good indication of how worn is the actuator. Typical minimum lifetime performance without abuse of the actuator is 10000 cycles in life test equals to 5,600,000 A x S.


### **Total number of over loads:**

Counts the number of times the lifting actuator has been overloaded.

### **Resetting of Service interval after service has been carried out:**

Resetting of service is done by pressing 2 buttons; Lifting arm UP, and lifting arm DOWN on the hand control at the same time for 5 seconds. After pressing the buttons for 5 seconds, you will receive an audio signal indicating that the timer has been reset. The timer will clear the display for service symbol, and start counting a new service period. This should only be done by an engineer or under guidance from a Prism Medical UK representative.

# Technical Specifications- Prism SA-200 Stand Aid



**Max.load and weight in kg and Dimensions in cm.**

Max. load	A		B	B'	C	D	F	F'	E	L		W1	W2	N.W.	Turning diameter
	Max.	Min.								Max.	Min.				
200	145.5	112	103.5	65	68~99	58~89	11	5	97.5	138	90	42.8	58.2	36.8	117

## DO'S & DON'TS

### DO

- ✓ Do charge the battery whenever possible. This will extend the battery life. A large number of cycles can be obtained from operating on the batteries, but battery lifetime is reduced with frequent discharging.
- ✓ Do inspect all cables particularly the mains power cable on the charger for any damage; replace where necessary.
- ✓ Do stow the handset and if fitted with charger, the mains power cable when transporting the hoist.
- ✓ Do clean the actuators, control box, charger, battery and handset at regular intervals to remove dust and dirt.

The exterior of the Stand Aid should only be cleaned, disinfected and sterilised using isopropyl alcohol. Damp a cloth with isopropyl alcohol and wipe down entire exterior of Stand Aid. No other chemicals and/or liquids should be used to clean, disinfect and sterilise this Stand Aid.

**⚠ Caution: Take great care to ensure that no liquids get inside the Stand Aid. This Stand Aid is not drip proof or water tight. Failure to protect the Stand Aid from liquids may result in damage to the Stand Aid and/or may cause personal injury.**

### DON'T

- ✗ Don't allow the batteries to fully discharge before connecting to the charger. The batteries are a lead-acid gel cell type that can be trickled charged continuously (batteries used for standby) and have a high current discharge capacity. Life is greatly reduced by deep or complete discharging of the batteries. Longer lifetime is obtained by maintaining fully charged batteries.
- ✗ Don't continue to operate the handset by repeatedly pressing the buttons if the lift function will not move or the actuator will not function. If this occurs then the actuator has either reached its end position, the load is too great or there is a malfunction.
- ✗ Don't exert excessive force on the handset cable as this may break off the wires inside the cable and prevent some or all of the operations.
- ✗ Don't continually operate the hoist functions. The system is not designed for continuous operation. Continuous Operation will cause batteries to be deeply discharged or to damage the actuator or control box by overheating.

# Fault Finding

Should problems arise with the use of the Stand Aid then review the following chart. Find the fault and complete the recommended solution. If the fault is not found and/or the solution does not correct the problem contact your local Prism Medical UK

PROBLEM	POSSIBLE CAUSE
Stand Aid cannot raise/lower, control box is not clicking.	<ul style="list-style-type: none"> <li>● Hand control is not properly plugged in.</li> <li>● Hand control is defective.</li> <li>● Control box is defective.</li> <li>● Emergency button is pushed in.</li> <li>● Mains cable is plugged into the control box.</li> </ul>
Stand Aid cannot raise/lower, control box is clicking.	<ul style="list-style-type: none"> <li>● Lifting actuator cord (1) is not fully plugged into slot (1) on the control box.</li> <li>● Lifting actuator is defective.</li> </ul>
Emergency Lower button on control box does not lower.	<ul style="list-style-type: none"> <li>● Lifting actuator cord (1) is not fully plugged into slot (1) on the control box.</li> <li>● Lifting actuator is defective.</li> <li>● Control box is defective.</li> </ul>
Stand Aid cannot spread or close legs.	<ul style="list-style-type: none"> <li>● Foot Pedal is defective.</li> </ul>
Cannot move Stand Aid.	<ul style="list-style-type: none"> <li>● Rear castor brake is engaged.</li> <li>● Castors are defective.</li> </ul>
Cannot tighten Mast into Base tube.	<ul style="list-style-type: none"> <li>● Mast locking assembly are defective, replace with new Mast locking assembly .</li> </ul>

# General Inspection and Maintenance

## A) Each Use - To be completed by User

- All functions on hand control are operational.

## B) Semi-Annually - To be completed at least every 6 months.

**Should any of these items fail the inspection do not use the Stand Aid. Contact Prism Medical UK or your local qualified service technician for service.**



- Complete the inspection as noted in the “Each Use” section above.

### **Also check the following:**

- Battery charger functional.
- Castors are clean and free of debris.
- All bolts and attachment points are secure (see next page for “Points of Attachment”)
- Visual inspection of all weld points and cast components.
- Check operation of the emergency stop button.

## C) Annually - To be completed at least every 12 months.



- Complete the inspection as noted in the “Semi-Annually” section above.

### **Also check the following:**


- Perform a working load test of one (1) lifting cycle with the maximum load. This is in accordance to the periodic inspection requirements of ISO10535.
- Check functional operation of the emergency lowering mechanism.
- All bolts in connection with the carry bar must be inspected, greased and tightened. All defective or worn components must be replaced.
- Connection bolts/bearings between arm and mast to be inspected. All defective or worn components must be replaced.
- Leg spreading function to be tested and all bolts to be inspected. All defective or worn components must be replaced.
- Check that main actuator is functional.



**This section to be only completed by a qualified service technician as authorised by Prism Medical UK . Use the provided service record history sheets to note findings during inspection.**

# SA200 Stand Aid- Points of Attachment

Inspect the bolts on a monthly basis in the areas shown.





# Service Record History - Initial Information

- Complete the following section on **Purchase and Service Information** as soon as this equipment is purchased.
- Use the service record history to record to any completed service and repairs.
- Ensure that the service record is signed and dated each time it is used.
- Be sure to have this piece of equipment serviced on a regular basis as described in the General Inspection and Maintenance Section.

## PURCHASE INFORMATION:

**Product Name:** SA-200 Stand Aid

**Model #:** \_\_\_\_\_

**Serial #:** \_\_\_\_\_

**Date of Purchase:** \_\_\_\_\_

**Purchased From:** \_\_\_\_\_

(local authorised Prism Medical UK Representative)

**Address:** \_\_\_\_\_

**City:** \_\_\_\_\_ **Postal Code:** \_\_\_\_\_

**Telephone No:** \_\_\_\_\_

**Comments:**

## SERVICE INFORMATION:

Contact the following company for service:

**Company:** \_\_\_\_\_

(local authorised Prism Medical UK Representative)

**Address:** \_\_\_\_\_

**City:** \_\_\_\_\_ **Postal Code:** \_\_\_\_\_

**Telephone No:** \_\_\_\_\_

**Comments:**

# Service Record History

Complete this section after each service, repair inspection and/or maintenance. Photocopy additional pages as required.

<b>Date:</b> _____		<b>Time:</b> _____	
<b>Service Type:</b> <input type="checkbox"/> Periodic Inspection <input type="checkbox"/> Monthly Inspection <input type="checkbox"/> 6 Month Inspection <input type="checkbox"/> Repair <input type="checkbox"/> Yearly Inspection <input type="checkbox"/> Other: _____			
<b>Completed By:</b> _____		_____	
Printed Name		Signature	
<b>Company:</b> _____			
<b>Remarks &amp; Action Taken:</b>			

<b>Date:</b> _____		<b>Time:</b> _____	
<b>Service Type:</b> <input type="checkbox"/> Periodic Inspection <input type="checkbox"/> Monthly Inspection <input type="checkbox"/> 6 Month Inspection <input type="checkbox"/> Repair <input type="checkbox"/> Yearly Inspection <input type="checkbox"/> Other: _____			
<b>Completed By:</b> _____		_____	
Printed Name		Signature	
<b>Company:</b> _____			
<b>Remarks &amp; Action Taken:</b>			

<b>Date:</b> _____		<b>Time:</b> _____	
<b>Service Type:</b> <input type="checkbox"/> Periodic Inspection <input type="checkbox"/> Monthly Inspection <input type="checkbox"/> 6 Month Inspection <input type="checkbox"/> Repair <input type="checkbox"/> Yearly Inspection <input type="checkbox"/> Other: _____			
<b>Completed By:</b> _____		_____	
Printed Name		Signature	
<b>Company:</b> _____			
<b>Remarks &amp; Action Taken:</b>			

<b>Date:</b> _____		<b>Time:</b> _____	
<b>Service Type:</b> <input type="checkbox"/> Periodic Inspection <input type="checkbox"/> Monthly Inspection <input type="checkbox"/> 6 Month Inspection <input type="checkbox"/> Repair <input type="checkbox"/> Yearly Inspection <input type="checkbox"/> Other: _____			
<b>Completed By:</b> _____		_____	
Printed Name		Signature	
<b>Company:</b> _____			
<b>Remarks &amp; Action Taken:</b>			

<b>Date:</b> _____		<b>Time:</b> _____	
<b>Service Type:</b> <input type="checkbox"/> Periodic Inspection <input type="checkbox"/> Monthly Inspection <input type="checkbox"/> 6 Month Inspection <input type="checkbox"/> Repair <input type="checkbox"/> Yearly Inspection <input type="checkbox"/> Other: _____			
<b>Completed By:</b> _____		_____	
Printed Name		Signature	
<b>Company:</b> _____			
<b>Remarks &amp; Action Taken:</b>			

<b>Date:</b> _____		<b>Time:</b> _____	
<b>Service Type:</b> <input type="checkbox"/> Periodic Inspection <input type="checkbox"/> Monthly Inspection <input type="checkbox"/> 6 Month Inspection <input type="checkbox"/> Repair <input type="checkbox"/> Yearly Inspection <input type="checkbox"/> Other: _____			
<b>Completed By:</b> _____		_____	
Printed Name		Signature	
<b>Company:</b> _____			
<b>Remarks &amp; Action Taken:</b>			

# Service Record History

Complete this section after each service, repair inspection and/or maintenance. Photocopy additional pages as required.

<b>Date:</b> _____		<b>Time:</b> _____	
<b>Service Type:</b> <input type="checkbox"/> Periodic Inspection <input type="checkbox"/> Monthly Inspection <input type="checkbox"/> 6 Month Inspection <input type="checkbox"/> Repair <input type="checkbox"/> Yearly Inspection <input type="checkbox"/> Other: _____			
<b>Completed By:</b> _____		_____	
Printed Name		Signature	
<b>Company:</b> _____			
<b>Remarks &amp; Action Taken:</b>			

<b>Date:</b> _____		<b>Time:</b> _____	
<b>Service Type:</b> <input type="checkbox"/> Periodic Inspection <input type="checkbox"/> Monthly Inspection <input type="checkbox"/> 6 Month Inspection <input type="checkbox"/> Repair <input type="checkbox"/> Yearly Inspection <input type="checkbox"/> Other: _____			
<b>Completed By:</b> _____		_____	
Printed Name		Signature	
<b>Company:</b> _____			
<b>Remarks &amp; Action Taken:</b>			

<b>Date:</b> _____		<b>Time:</b> _____	
<b>Service Type:</b> <input type="checkbox"/> Periodic Inspection <input type="checkbox"/> Monthly Inspection <input type="checkbox"/> 6 Month Inspection <input type="checkbox"/> Repair <input type="checkbox"/> Yearly Inspection <input type="checkbox"/> Other: _____			
<b>Completed By:</b> _____		_____	
Printed Name		Signature	
<b>Company:</b> _____			
<b>Remarks &amp; Action Taken:</b>			

<b>Date:</b> _____		<b>Time:</b> _____	
<b>Service Type:</b> <input type="checkbox"/> Periodic Inspection <input type="checkbox"/> Monthly Inspection <input type="checkbox"/> 6 Month Inspection <input type="checkbox"/> Repair <input type="checkbox"/> Yearly Inspection <input type="checkbox"/> Other: _____			
<b>Completed By:</b> _____		_____	
Printed Name		Signature	
<b>Company:</b> _____			
<b>Remarks &amp; Action Taken:</b>			

<b>Date:</b> _____		<b>Time:</b> _____	
<b>Service Type:</b> <input type="checkbox"/> Periodic Inspection <input type="checkbox"/> Monthly Inspection <input type="checkbox"/> 6 Month Inspection <input type="checkbox"/> Repair <input type="checkbox"/> Yearly Inspection <input type="checkbox"/> Other: _____			
<b>Completed By:</b> _____		_____	
Printed Name		Signature	
<b>Company:</b> _____			
<b>Remarks &amp; Action Taken:</b>			

<b>Date:</b> _____		<b>Time:</b> _____	
<b>Service Type:</b> <input type="checkbox"/> Periodic Inspection <input type="checkbox"/> Monthly Inspection <input type="checkbox"/> 6 Month Inspection <input type="checkbox"/> Repair <input type="checkbox"/> Yearly Inspection <input type="checkbox"/> Other: _____			
<b>Completed By:</b> _____		_____	
Printed Name		Signature	
<b>Company:</b> _____			
<b>Remarks &amp; Action Taken:</b>			

Copyright Prism Medical UK 2014

# Service Record History

Complete this section after each service, repair inspection and/or maintenance. Photocopy additional pages as required.

<b>Date:</b> _____		<b>Time:</b> _____	
<b>Service Type:</b> <input type="checkbox"/> Periodic Inspection <input type="checkbox"/> Monthly Inspection <input type="checkbox"/> 6 Month Inspection <input type="checkbox"/> Repair <input type="checkbox"/> Yearly Inspection <input type="checkbox"/> Other: _____			
<b>Completed By:</b> _____		_____	
Printed Name		Signature	
<b>Company:</b> _____			
<b>Remarks &amp; Action Taken:</b>			

<b>Date:</b> _____		<b>Time:</b> _____	
<b>Service Type:</b> <input type="checkbox"/> Periodic Inspection <input type="checkbox"/> Monthly Inspection <input type="checkbox"/> 6 Month Inspection <input type="checkbox"/> Repair <input type="checkbox"/> Yearly Inspection <input type="checkbox"/> Other: _____			
<b>Completed By:</b> _____		_____	
Printed Name		Signature	
<b>Company:</b> _____			
<b>Remarks &amp; Action Taken:</b>			

<b>Date:</b> _____		<b>Time:</b> _____	
<b>Service Type:</b> <input type="checkbox"/> Periodic Inspection <input type="checkbox"/> Monthly Inspection <input type="checkbox"/> 6 Month Inspection <input type="checkbox"/> Repair <input type="checkbox"/> Yearly Inspection <input type="checkbox"/> Other: _____			
<b>Completed By:</b> _____		_____	
Printed Name		Signature	
<b>Company:</b> _____			
<b>Remarks &amp; Action Taken:</b>			

<b>Date:</b> _____		<b>Time:</b> _____	
<b>Service Type:</b> <input type="checkbox"/> Periodic Inspection <input type="checkbox"/> Monthly Inspection <input type="checkbox"/> 6 Month Inspection <input type="checkbox"/> Repair <input type="checkbox"/> Yearly Inspection <input type="checkbox"/> Other: _____			
<b>Completed By:</b> _____		_____	
Printed Name		Signature	
<b>Company:</b> _____			
<b>Remarks &amp; Action Taken:</b>			

<b>Date:</b> _____		<b>Time:</b> _____	
<b>Service Type:</b> <input type="checkbox"/> Periodic Inspection <input type="checkbox"/> Monthly Inspection <input type="checkbox"/> 6 Month Inspection <input type="checkbox"/> Repair <input type="checkbox"/> Yearly Inspection <input type="checkbox"/> Other: _____			
<b>Completed By:</b> _____		_____	
Printed Name		Signature	
<b>Company:</b> _____			
<b>Remarks &amp; Action Taken:</b>			

<b>Date:</b> _____		<b>Time:</b> _____	
<b>Service Type:</b> <input type="checkbox"/> Periodic Inspection <input type="checkbox"/> Monthly Inspection <input type="checkbox"/> 6 Month Inspection <input type="checkbox"/> Repair <input type="checkbox"/> Yearly Inspection <input type="checkbox"/> Other: _____			
<b>Completed By:</b> _____		_____	
Printed Name		Signature	
<b>Company:</b> _____			
<b>Remarks &amp; Action Taken:</b>			

# Service Record History

Complete this section after each service, repair inspection and/or maintenance. Photocopy additional pages as required.

<b>Date:</b> _____		<b>Time:</b> _____	
<b>Service Type:</b> <input type="checkbox"/> Periodic Inspection <input type="checkbox"/> Monthly Inspection <input type="checkbox"/> 6 Month Inspection <input type="checkbox"/> Repair <input type="checkbox"/> Yearly Inspection <input type="checkbox"/> Other: _____			
<b>Completed By:</b> _____		_____	
Printed Name		Signature	
<b>Company:</b> _____			
<b>Remarks &amp; Action Taken:</b>			

<b>Date:</b> _____		<b>Time:</b> _____	
<b>Service Type:</b> <input type="checkbox"/> Periodic Inspection <input type="checkbox"/> Monthly Inspection <input type="checkbox"/> 6 Month Inspection <input type="checkbox"/> Repair <input type="checkbox"/> Yearly Inspection <input type="checkbox"/> Other: _____			
<b>Completed By:</b> _____		_____	
Printed Name		Signature	
<b>Company:</b> _____			
<b>Remarks &amp; Action Taken:</b>			

<b>Date:</b> _____		<b>Time:</b> _____	
<b>Service Type:</b> <input type="checkbox"/> Periodic Inspection <input type="checkbox"/> Monthly Inspection <input type="checkbox"/> 6 Month Inspection <input type="checkbox"/> Repair <input type="checkbox"/> Yearly Inspection <input type="checkbox"/> Other: _____			
<b>Completed By:</b> _____		_____	
Printed Name		Signature	
<b>Company:</b> _____			
<b>Remarks &amp; Action Taken:</b>			

<b>Date:</b> _____		<b>Time:</b> _____	
<b>Service Type:</b> <input type="checkbox"/> Periodic Inspection <input type="checkbox"/> Monthly Inspection <input type="checkbox"/> 6 Month Inspection <input type="checkbox"/> Repair <input type="checkbox"/> Yearly Inspection <input type="checkbox"/> Other: _____			
<b>Completed By:</b> _____		_____	
Printed Name		Signature	
<b>Company:</b> _____			
<b>Remarks &amp; Action Taken:</b>			

<b>Date:</b> _____		<b>Time:</b> _____	
<b>Service Type:</b> <input type="checkbox"/> Periodic Inspection <input type="checkbox"/> Monthly Inspection <input type="checkbox"/> 6 Month Inspection <input type="checkbox"/> Repair <input type="checkbox"/> Yearly Inspection <input type="checkbox"/> Other: _____			
<b>Completed By:</b> _____		_____	
Printed Name		Signature	
<b>Company:</b> _____			
<b>Remarks &amp; Action Taken:</b>			

<b>Date:</b> _____		<b>Time:</b> _____	
<b>Service Type:</b> <input type="checkbox"/> Periodic Inspection <input type="checkbox"/> Monthly Inspection <input type="checkbox"/> 6 Month Inspection <input type="checkbox"/> Repair <input type="checkbox"/> Yearly Inspection <input type="checkbox"/> Other: _____			
<b>Completed By:</b> _____		_____	
Printed Name		Signature	
<b>Company:</b> _____			
<b>Remarks &amp; Action Taken:</b>			

Copyright Prism Medical UK 2014

# Service Record History

Complete this section after each service, repair inspection and/or maintenance. Photocopy additional pages as required.

<b>Date:</b> _____		<b>Time:</b> _____	
<b>Service Type:</b> <input type="checkbox"/> Periodic Inspection <input type="checkbox"/> Monthly Inspection <input type="checkbox"/> 6 Month Inspection <input type="checkbox"/> Repair <input type="checkbox"/> Yearly Inspection <input type="checkbox"/> Other: _____			
<b>Completed By:</b> _____		_____	
Printed Name		Signature	
<b>Company:</b> _____			
<b>Remarks &amp; Action Taken:</b>			

<b>Date:</b> _____		<b>Time:</b> _____	
<b>Service Type:</b> <input type="checkbox"/> Periodic Inspection <input type="checkbox"/> Monthly Inspection <input type="checkbox"/> 6 Month Inspection <input type="checkbox"/> Repair <input type="checkbox"/> Yearly Inspection <input type="checkbox"/> Other: _____			
<b>Completed By:</b> _____		_____	
Printed Name		Signature	
<b>Company:</b> _____			
<b>Remarks &amp; Action Taken:</b>			

<b>Date:</b> _____		<b>Time:</b> _____	
<b>Service Type:</b> <input type="checkbox"/> Periodic Inspection <input type="checkbox"/> Monthly Inspection <input type="checkbox"/> 6 Month Inspection <input type="checkbox"/> Repair <input type="checkbox"/> Yearly Inspection <input type="checkbox"/> Other: _____			
<b>Completed By:</b> _____		_____	
Printed Name		Signature	
<b>Company:</b> _____			
<b>Remarks &amp; Action Taken:</b>			

<b>Date:</b> _____		<b>Time:</b> _____	
<b>Service Type:</b> <input type="checkbox"/> Periodic Inspection <input type="checkbox"/> Monthly Inspection <input type="checkbox"/> 6 Month Inspection <input type="checkbox"/> Repair <input type="checkbox"/> Yearly Inspection <input type="checkbox"/> Other: _____			
<b>Completed By:</b> _____		_____	
Printed Name		Signature	
<b>Company:</b> _____			
<b>Remarks &amp; Action Taken:</b>			

<b>Date:</b> _____		<b>Time:</b> _____	
<b>Service Type:</b> <input type="checkbox"/> Periodic Inspection <input type="checkbox"/> Monthly Inspection <input type="checkbox"/> 6 Month Inspection <input type="checkbox"/> Repair <input type="checkbox"/> Yearly Inspection <input type="checkbox"/> Other: _____			
<b>Completed By:</b> _____		_____	
Printed Name		Signature	
<b>Company:</b> _____			
<b>Remarks &amp; Action Taken:</b>			

<b>Date:</b> _____		<b>Time:</b> _____	
<b>Service Type:</b> <input type="checkbox"/> Periodic Inspection <input type="checkbox"/> Monthly Inspection <input type="checkbox"/> 6 Month Inspection <input type="checkbox"/> Repair <input type="checkbox"/> Yearly Inspection <input type="checkbox"/> Other: _____			
<b>Completed By:</b> _____		_____	
Printed Name		Signature	
<b>Company:</b> _____			
<b>Remarks &amp; Action Taken:</b>			

# NOTES

# Warranty

This Warranty does not affect or in any way limit your Statutory Rights

- 1) Prism Medical UK guarantees all equipment supplied as new, against failure within the period of 1 year from date of purchase by virtue of defects in material or workmanship.
- 2) This warranty does not apply to failure attributable to normal wear and tear, damage by natural forces, user neglect or misuse or to deliberate destruction, or to batteries more than 90 days after original purchase.
- 3) This equipment warranty shall be void if the equipment is not serviced by Prism Medical UK or its authorised service agents in accordance with the manufacturer's recommendations or if any unauthorised person carries out works on the equipment.
- 4) The liability of Prism Medical UK under the terms of this guarantee shall be limited to the replacement of defective part (s) to the sales distributor, dealer, agent, person or entity which purchased the equipment from Prism Medical UK. In no event shall Prism Medical UK incur liability for any consequential or unforeseeable losses.

---

If you have any questions about the manufacture or operation of this equipment, please contact Prism Medical UK, or your local authorised dealer.

**This document conforms to EN ISO 10535 requirements**

**PRISM MEDICAL UK**

Prism Medical UK, Unit 4 Jubilee Business Park, Jubilee Way,  
Grange Moor, Wakefield, WF4 4TD

Tel 0845 409 2000

[www.prismmedical.co.uk](http://www.prismmedical.co.uk)

[info@prismmedical.co.uk](mailto:info@prismmedical.co.uk)

