

# Shower Chair Technical

Below are the technical specifications of the most popular Freeway Shower Chair models. Please consult our website, [www.prismmedical.co.uk](http://www.prismmedical.co.uk) for the full range of dimensions.

- 1 = Overall maximum height.
- 2 = Overall maximum width.
- 3 = Length including footrests (if fitted) at lowest position.
- 4 = Length including footrests (if fitted) at highest position.
- 5 = Width of standard seat.
- 6 = Height to top of seat.
- 7 = Height clearance.
- 8 = Seat depth (standard seat).
- 9 = Top of armrest to top of seat.
- 10 = Width clearance.
- 11 = Top of backrest to seat.
- 12 = Distance between armrests at rear.
- 13 = Distance between armrests at front.

T30													
	1	2	3	4	5	6	7	8	9	10	11	12	13
T30A (SC only)	920	640	570	570	540	460	345	400	180	485	450	470	530
T30B (SC only)	920	590	570	570	490	460	345	400	180	435	450	420	480
T30C (SC only)	890	540	570	570	440	460	345	400	180	385	450	370	430

T40													
	1	2	3	4	5	6	7	8	9	10	11	12	13
T40A	980	640	1050	960	540	560	445	400	180	485	450	470	530
T40B	980	590	1050	960	490	560	445	400	180	435	450	420	480
T40C	980	540	1050	960	440	560	445	400	180	385	450	370	430
T40A AUTO	1000	590	1050	960	540	560	480	400	180	485	450	460	530
T40B AUTO	1000	540	1050	960	490	560	480	400	180	435	450	410	480
T40C AUTO	1000	540	1050	960	440	560	480	400	180	385	450	410	430
T40A COM	980	690	1050	960	540	560	470	400	180	485	450	470	530
T40B COM	980	590	1050	960	490	560	470	400	180	435	450	420	480
T40C COM	980	540	1050	960	440	560	470	400	180	385	450	370	430

T50													
	1	2	3	4	5	6	7	8	9	10	11	12	13
T50A	980	730	940	900	540	560	445	400	180	485	450	470	530
T50B	980	680	940	900	490	560	445	400	180	435	450	420	480
T50C	980	630	940	900	440	560	445	400	180	385	450	370	430

T60													
	1	2	3	4	5	6	7	8	9	10	11	12	13
T60A	980	710	1050	970	540	550	435	400	180	485	450	470	530
T60B	980	660	1050	970	490	550	435	400	180	435	450	420	480
T60C	980	610	1050	970	440	550	435	400	180	385	450	370	430

T70													
	1	2	3	4	5	6	7	8	9	10	11	12	13
T70A	980	750	1150	1100	540	560	445	400	180	485	450	470	530
T70B	980	700	1150	1100	490	560	445	400	180	435	450	420	480
T70C	980	650	1150	1100	440	560	445	400	180	385	450	370	430
T70A AUTO	1000	790	1140	1060	540	560	480	400	180	485	450	460	530
T70B AUTO	1000	740	1140	1060	490	560	480	400	180	435	450	410	480
T70C AUTO	1000	690	1140	1060	440	560	480	400	180	385	450	410	430
T70A COM	980	775	1110	1020	540	560	470	400	180	485	450	470	530
T70B COM	980	725	1110	1020	490	560	470	400	180	435	450	420	480
T70C COM	980	675	1110	1020	440	560	470	400	180	385	450	370	430

T90 Paediatric													
	1	2	3	4	5	6	7	8	9	10	11	12	13
T90 Assistant Prop.	920	620	730	N/A	360	460	N/A	330	130	N/A	440	300	340
T90 Self Propelled	920	620	880	N/A	360	460	N/A	330	130	N/A	440	300	340
T90WC Assistant Pr.	980	470	680	N/A	360	560	470	330	130	430	440	320	320
T90WC Self Prop.	980	680	880	N/A	360	560	470	330	130	430	440	320	320

T100 Bariatric													
	1	2	3	4	5	6	7	8	9	10	11	12	13
T100	990	900	1220	1170	790	570	450	550	170	730	430	720	760
T100FR	990	900	1220	1170	790	570	450	550	170	730	430	720	760

T80 Recliner															
	1	2	3	4A	4B	4C	5	6	7	8	9	10	11	12	13
T80A	1125	1440	640	1040	1300	1500	540	570	445	400	190	480	535	470	530
T80B	1125	1440	590	1040	1300	1500	490	570	445	400	190	430	535	420	480
T80C	1125	1440	540	1040	1300	1500	440	570	445	400	190	380	535	370	430
T80A COM	1125	1440	730	1040	1300	1500	540	570	470	400	190	480	535	460	530
T80B COM	1125	1440	680	1040	1300	1500	490	570	470	400	190	430	535	410	480
T80C COM	1125	1440	630	1040	1300	1500	440	570	470	400	190	380	535	410	430
T80A AUTO	1125	1440	730	1040	1300	1500	540	570	470	400	190	480	535	470	530
T80B AUTO	1125	1440	680	1040	1300	1500	490	570	470	400	190	430	535	420	480
T80C AUTO	1125	1440	630	1040	1300	1500	440	570	470	400	190	430	535	370	430
T80 Paediatric	1050	1345	590	760	1050	1170	360	570	470	325	145	430	530	360	360

**T80 Specification Key**

- 1 = Height.
- 2 = Height with detachable headrest.
- 3 = Width.
- 4A = Length in upright position.

- 4B = Length with max. recline including footrests at highest position.
- 4C = Length max. recline and detachable headrest fitted.
- 5 = Seat width.
- 6 = Height to top of seat.
- 7 = Height clearance.
- 8 = Seat depth (standard seat).

- 9 = Top of armrest to top of seat.
- 10 = Width clearance.
- 11 = Top of backrest to top of seat (without detachable headrest fitted - with headrest add 300mm).
- 12 = Distance between armrests at rear.
- 13 = Distance between armrests at front.

Contact Us: To request a quote or assessment, please call us on **0845 409 2000** or visit [www.prismmedical.co.uk](http://www.prismmedical.co.uk)