

Etac 4Way Glide



PRISM MEDICAL UK

Moving, Handling and Bathing Solutions™

ETAC 4Way Glide

The Etac 4Way Glide allows easy low friction movement to provide a permanent 'on the bed' solution for patients who require frequent repositioning.

The complete 4Way Glide system uses two components, a nylon undersheet and a top quilt with either breathable PU or cotton surface.



Key Features

- Reduces the incidence of pressure ulcers
- Demands minimal effort from carers
- Ideal for use by a single carer
- Works flawlessly with four section profiling beds
- Can be used on a pressure care mattress
- May be used with a hoist to assist turning
- Top quilt machine washable to 80°C / Nylon sheet can be washed up to 95°C

Options Available

	ORDER CODES		
Product	Standard	Bariatric	Paediatric
4 Way Complete (cotton)	4WAYCOMPCOT	BAR4WAYCOMPCO	PED4WAYCOMPCOT
4 Way Complete (PU)	4WAYCOMPPU	BAR4WAYCOMPPU	PED4WAYCOMPPU
Nylon Undersheet	4WAYNYCOV	BAR4WAYNYCOV	PED4WAYNYCOV
Cotton Top Quilt	4WAYQUILT	BAR4WAYQUILT	PED4WAYQUILT
PU Top Quilt	4WAYQUILTPU	BAR4WAYQUILTPU	PED4WAYQUILTPU

To arrange a demonstration of the 4Way Glide please contact us on 0845 409 2000 or visit www.prismmedical.co.uk