

A new cost effective range to cater for a variety of transfer and hoisting requirements

Mackworth SA160E

The **Mackworth SA160E** offers a highly economical stand aid solution.

The **SA160E** stand aid is produced to the same cost effective manufacturing process of the Mackworth hoist range, providing a robust unit with **no compromise on quality**.

Offers true value-for-money.

Key Benefits and Features:

- A cost effective solution for the budget conscious
- Accommodates the support and comfort needs of the client
- Maximum weight capacity – 160kg
- Manual leg opening & closing
- Lifting range of 770 – 1,675mm
- Stand aid foot print (L x W) – 1,250 x 660mm
- Turning radius – 1,305mm
- One year warranty

Additional Technical Information

Legs Open – External Width – 1,190mm

Legs Open – Internal Width – 1,030mm

Legs Closed – External Width – 720mm

Legs Closed – Internal Width – 560mm

Overall Height of Legs (Top of Legs to Floor) – 112mm

Ground Clearance (Underside of Leg to Floor) – 62mm

Front Twin Castors (Diameter) – 100mm

Rear Twin Castors (Diameter) – 100mm

Weight Assembled – 46kg

Mast & Boom (including all fittings) – 23.1kg

Base Assembly – 22.8kg



**INTRODUCTORY PRICE OFFER
AVAILABLE NOW.
PLEASE SPEAK TO OUR EXPERT
ADVISERS.**

01924 840 100

www.mackworth-healthcare.com



Manufactured in the UK

Mackworth

Part of **PRISM MEDICAL UK**



A new cost effective range to cater for a variety of transfer and hoisting requirements

Mackworth TA180

The **Mackworth TA180** is a transfer support unit that encourages more mobile clients to stand un-assisted.

The versatility of the **Mackworth TA180** allows for ease of use in a variety of care and health environments.

Key Benefits and Features:

- Provides easy transportation to the bathroom and toilet
- Offers a suitable alternative to a wheelchair
- Allows for greater interaction between user and carer
- Padded knee support for extra comfort
- Maximum weight capacity – 180kg
- Stand aid foot print (L x W) – 1,500 x 650mm
- Manual leg open & closing
- Turning radius – 1,125mm
- One year warranty

Additional Technical Information

Legs Open – External Width – 1,600mm

Legs Open – Internal Width – 910mm

Legs Closed – External Width – 665mm

Legs Closed – Internal Width – 515mm

Overall Height of Legs (Top of Legs to Floor) – 113mm

Ground Clearance (Underside of Leg to Floor) – 64mm

Front and Rear Castors (Diameter) – 100mm

Rear Twin Castors (Diameter) – 100mm

Weight Assembled – 37kg



INTRODUCTORY PRICE OFFER
AVAILABLE NOW.
PLEASE SPEAK TO OUR EXPERT
ADVISERS.

01924 840 100

www.mackworth-healthcare.com



Manufactured in the UK

Mackworth

Part of PRISM MEDICAL UK

